

Antibodies Product Datasheet

Product Name:	eEF2 mouse mAb
Catalog Number:	N0860
Reactivity:	Human
Applications:	WB, ICC
Dilutions:	wb 1:5000 icc 1:200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human eEF2 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	eef2
Protein Name:	
Other Name:	EEF2; Eef2; EF-2; EF2; EF2_HUMAN; Elongation factor 2; Eukaryotic translation elongation factor 2; Polypeptidyl tRNA translocase; SCA26.
Human Swiss Prot No.:	P13639
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of eEF2 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year