

Antibodies Product Datasheet

Product Name:	Hexokinase II mouse mAb
Catalog Number:	N0863
Reactivity:	Rat, Mouse, Monkey
Applications:	WB
Dilutions:	wb 1:500
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant human Hexokinase II protein.
Concentration:	1 mg/ml
Gene Name:	hk2
Protein Name:	
Other Name:	DKFZp686M1669; Hexokinase 2; Hexokinase 2 muscle; Hexokinase type II; Hexokinase-2; HK 2; HK II; HK2; HKII; HxK 2; HxK2; HXK2_HUMAN; Muscle form hexokinase.
Human Swiss Prot No.:	P52789
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Hexokinase II and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year