

Antibodies Product Datasheet

Product Name:	CHRD1 mouse mAb
Catalog Number:	N0868
Reactivity:	Human
Applications:	WB,IHC-p,ICC
Dilutions:	wb 1:1000 IHC-p 1:300 icc 1:300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human CHRD1 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	chrd1
Protein Name:	
Other Name:	CHL;chordin like 1;Chordin-like protein 1;CHRD1;CRDL1_HUMAN;dA141H5.1;neuralin 1; Neuralin-1;Neurogenesin-1;NRLN1;Ventroptin;VOPT.
Human Swiss Prot No.:	Q9BU40
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CHRD1 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year