

Antibodies Product Datasheet

Product Name:	RanBP9 mouse mAb
Catalog Number:	N0869
Reactivity:	Human
Applications:	WB, ICC, IP
Dilutions:	wb 1:1000 icc 1:300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human RanBP9 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	ranbp9
Protein Name:	
Other Name:	B cell antigen receptor Ig beta associated protein 1;BPM 90;BPM L;BPM-L;BPM90;BPML;IBAP 1;Imp 9;Importin 9;Novel centrosomal protein RanBPM;RAN binding protein 9;Ran binding protein centrosomal;Ran Binding Protein in the Microtubule organizing center;Ran
Human Swiss Prot No.:	Q96S59
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of RanBP9 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year