

Antibodies Product Datasheet

Product Name:	hnRNP U mouse mAb
Catalog Number:	N0870
Reactivity:	Human
Applications:	WB
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human hnRNP U protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	hnrnpu
Protein Name:	
Other Name:	Heterogeneous nuclear ribonucleoprotein U;hnRNP U protein;HNRNPU;hnRNPU protein;HNRPU; HNRPU_HUMAN;p120;p120 nuclear protein;SAF A;SAF-A;SAFA;Scaffold attachment factor A;U21.1.
Human Swiss Prot No.:	Q00839
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of hnRNP U and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year