

Antibodies Product Datasheet

Product Name:	β-Arrestin 1 mouse mAb
Catalog Number:	N0871
Reactivity:	Mouse,Rat,Human
Applications:	WB
Dilutions:	wb 1:500
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human beta Arrestin 1 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	arrb1
Protein Name:	
Other Name:	ARB 1;ARB 1;ARB1;ARB1;ARR 1;ARR 1;ARR1;ARR1;ARRB 1;ARRB1;ARRB1;ARRB1_HUMAN;Arrestin 2; Arrestin beta 1;Arrestin beta-1;Beta-arrestin-1.
Human Swiss Prot No.:	P49407
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of beta Arrestin 1 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year