

Antibodies Product Datasheet

Product Name:	Fyn mouse mAb
Catalog Number:	N0872
Reactivity:	Human, Monkey, Rat, Mouse
Applications:	WB, ICC
Dilutions:	wb 1:500 icc 1:50
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human Fyn protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	fyn
Protein Name:	
Other Name:	C syn protooncogene; Fyn; FYN oncogene related to SRC FGR YES; FYN_HUMAN; MGC45350; OKT3 induced calcium influx regulator; P59 FYN; p59-Fyn; Protein tyrosine kinase fyn; Proto oncogene tyrosine protein kinase fyn; Proto-oncogene c-Fyn; Proto-oncogene Syn; Protooncoge
Human Swiss Prot No.:	P06241
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Fyn and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year