

Antibodies Product Datasheet

Product Name:	c-JUN mouse mAb
Catalog Number:	N0874
Reactivity:	
Applications:	WB
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human c-JUN protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	jun
Protein Name:	
Other Name:	Activator Protein 1;AP 1;AP1;cJun;Enhancer Binding Protein AP1;Jun Activation Domain Binding Protein;JUN;Jun oncogene;JUN protein;Jun proto oncogene;JUN_HUMAN;JUNC;Oncogene JUN;p39;Proto oncogene c jun;Proto oncogene cJun;Proto-oncogene c-jun;Transcriptio
Human Swiss Prot No.:	P05412
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Transfected
Storage and Stability:	-20°C/1 year