

Antibodies Product Datasheet

Product Name:	DNA-PKcs mouse mAb
Catalog Number:	N0877
Reactivity:	Human
Applications:	WB,IP,ICC,IHC
Dilutions:	wb 1:1000 icc 1:100
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human DNA-PKcs protein fragments expressed in E.coli
Concentration:	1 mg/ml
Gene Name:	prkdc
Protein Name:	
Other Name:	DNA dependent protein kinase catalytic subunit;DNA PKcs;DNA-dependent protein kinase catalytic subunit;DNA-PK catalytic subunit;DNAPK;DNPK1;hyper radiosensitivity of murine scid mutation, complementing 1;HYRC;HYRC1;p350;p460;PRKDC;PRKDC_HUMAN;protein kinase
Human Swiss Prot No.:	P78527
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DNA-PKcs and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year