

Antibodies Product Datasheet

Product Name:	Aldehyde dehydrogenase 10 mouse mAb
Catalog Number:	N0879
Reactivity:	Human
Applications:	WB
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human Aldehyde dehydrogenase 10 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	aldh3a2
Protein Name:	
Other Name:	Aldehyde dehydrogenase;Aldehyde dehydrogenase 10;Aldehyde dehydrogenase 3 family member A2;ALDH10;ALDH3A2;DKFZp686E23276;FALDH;Fatty aldehyde dehydrogenase;FLJ20851;Microsomal aldehyde dehydrogenase;OTTHUMP00000065799;OTTHUMP00000065800;OTT HUMP00000065801
Human Swiss Prot No.:	P51648
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Aldehyde dehydrogenase 10 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year

N0879 | Aldehyde dehydrogenase 10 mouse mAb

Tech@assaybiotech.com | Order@assaybiotech.com

Research use only not intended for therapeutic/diagnostic procedures

Last updated December 2022