

Antibodies Product Datasheet

Product Name:	Jak1 mouse mAb
Catalog Number:	N0881
Reactivity:	Human,Rat
Applications:	WB,IF/ICC,IP
Dilutions:	wb 1:1000 icc 1:200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human Jak1 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	jak1
Protein Name:	
Other Name:	JAK 1; JAK 1A; JAK 1B; JAK-1; JAK1; JAK1_HUMAN; JAK1A; JAK1B; Janus kinase 1 (a protein tyrosine kinase); Janus kinase 1; JTK3; Tyrosine protein kinase JAK 1; Tyrosine protein kinase JAK1; Tyrosine-protein kinase JAK1.
Human Swiss Prot No.:	P23458
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Jak1 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year