

Antibodies Product Datasheet

Product Name:	LRRFIP1 mouse mAb
Catalog Number:	N0885
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human LRRFIP1 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	lrrfip1
Protein Name:	
Other Name:	FLAP1; FLIAP 1; GC binding factor 2; GCF2; HUF11; Leucine rich repeat (in FLII) interacting protein 1; Leucine rich repeat flightless interacting protein 1; LRR FLII interacting protein 1; LRRFIP 1; Lrrfip1; MGC10947; MGC119738; MGC119739; NEDD8 conjugating enzyme; TA
Human Swiss Prot No.:	Q32MZ4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of LRRFIP1 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year