

Antibodies Product Datasheet

Product Name:	Calumenin mouse mAb
Catalog Number:	N0890
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human Calumenin protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	calu
Protein Name:	
Other Name:	CALU; CALU_HUMAN; Calumenin 1; Calumenin 2; Calumenin; Crocalbin; crocalbin-like protein; FLJ90608; IEF SSP 9302; multiple EF-hand protein; OTTHUMP00000212712; OTTHUMP00000212713; OTTHUMP00000212714; OTTHUMP00000212715.
Human Swiss Prot No.:	O43852
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Calumenin and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year