

Antibodies Product Datasheet

Product Name:	Beclin 1 mouse mAb
Catalog Number:	N0891
Reactivity:	Human, Mouse, Rat, Monkey, Hamster
Applications:	WB, IP
Dilutions:	wb 1:1000 IP:1:50-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human Beclin 1 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	becn1
Protein Name:	
Other Name:	APG6; ATG 6; ATG6; ATG6 autophagy related 6 homolog; Bcl-2-interacting protein beclin; Beclin 1 (coiled coil moesin like BCL2 interacting protein); Beclin 1 autophagy related; Beclin-1; Beclin1; BECN1; BECN1; BECN1_HUMAN; Coiled coil myosin like BCL2 interacting pro
Human Swiss Prot No.:	Q14457
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Beclin 1 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year