

Antibodies Product Datasheet

Product Name:	LCK mouse mAb
Catalog Number:	N0892
Reactivity:	Human
Applications:	WB
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human LCK protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	lck
Protein Name:	
Other Name:	LCK;Lck p56;LCK_HUMAN;Leukocyte C-terminal Src kinase;LSK;Lymphocyte cell specific protein tyrosine kinase;Lymphocyte cell-specific protein-tyrosine kinase;Lymphocyte Specific Protein Tyrosine Kinase;Membrane associated protein tyrosine kinase;Oncogene lc
Human Swiss Prot No.:	P06239
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of LCK and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year