

Antibodies Product Datasheet

Product Name:	Pyruvate Dehydrogenase E2 mouse mAb
Catalog Number:	N0893
Reactivity:	Human, Mouse
Applications:	WB, ICC, IP
Dilutions:	wb 1:1000 icc 1:300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human Pyruvate Dehydrogenase E2 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	dlat
Protein Name:	
Other Name:	70 kDa mitochondrial autoantigen of primary biliary cirrhosis; anti DLAT; Dihydrolipoamide acetyltransferase component of pyruvate dehydrogenase complex; Dihydrolipoamide; Dihydrolipoamide S Acetyltransferase; Dihydrolipoamide S-acetyltransferase (E2 compone
Human Swiss Prot No.:	P10515
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Pyruvate Dehydrogenase E2 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year