

Antibodies Product Datasheet

Product Name:	AKR1B1 mouse mAb
Catalog Number:	N0896
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human AKR1B1 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	akr1b1
Protein Name:	
Other Name:	ADR; AKR1B 1; Akr1b1; Aldehyde reductase 1; Aldehyde reductase; Aldo keto reductase family 1, member B1; Aldo-keto reductase family 1 member B1; aldo-keto reductase family 1, member B1 (aldose reductase); Aldose reductase; aldr 1; ALDR_HUMAN; ALDR1; ALR2; AR; Lii5 2 C
Human Swiss Prot No.:	P15121
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of AKR1B1 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year