

Antibodies Product Datasheet

Product Name:	TXNRD2 mouse mAb
Catalog Number:	N0897
Reactivity:	Human
Applications:	WB
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human TXNRD2 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	txnrd2
Protein Name:	
Other Name:	mitochondrial;selenoprotein Z;SELZ;Thioredoxin reductase 2;Thioredoxin reductase 2 mitochondrial;thioredoxin reductase 3;thioredoxin reductase beta;Thioredoxin reductase TR3;TR 3;TR;TR beta;TR-beta;TRXR 2;TRXR2;TRXR2_HUMAN;TXNRD 2;Txnrd2.
Human Swiss Prot No.:	Q9NNW7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of TXNRD2 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year