

Antibodies Product Datasheet

Product Name:	Akt1 mouse mAb
Catalog Number:	N0899
Reactivity:	Mouse, Transfected
Applications:	WB
Dilutions:	wb 1:2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human AKT1 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	akt1
Protein Name:	
Other Name:	AKT 1; AKT; AKT1; AKT1_HUMAN; C AKT; MGC99656; PKB; PKB-ALPHA; PRKBA; Protein Kinase B Alpha; Protein kinase B; Proto-oncogene c-Akt; RAC Alpha; RAC; RAC PK Alpha; RAC Serine/Threonine Protein Kinase; RAC-alpha serine/threonine-protein kinase; RAC-PK-alpha; vAKT Murine Thy
Human Swiss Prot No.:	P31749
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of AKT1 and does not cross-react with AKT2 and other proteins.
Storage and Stability:	-20°C/1 year