

Antibodies Product Datasheet

Product Name:	HDAC3 mouse mAb
Catalog Number:	N0900
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human HDAC3 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	hdac3
Protein Name:	
Other Name:	HD 3; HD3; HDAC 3; HDAC3; HDAC3_HUMAN; Histone deacetylase 3 (HD3) (RPD3-2); histone deacetylase 3; RPD 3; RPD3 2; RPD3; RPD3-2; SMAP 45; SMAP45.
Human Swiss Prot No.:	O15379
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of HDAC3 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year