

Antibodies Product Datasheet

Product Name:	Musashi-1 mouse mAb
Catalog Number:	N0901
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human Musashi-1 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	msi1
Protein Name:	
Other Name:	MSI1; Musashi-1; Msi 1; Msi1; MSI1H_HUMAN; Musashi homolog 1; Musashi-1; Musashi1; RNA binding protein Musashi homolog 1; RNA-binding protein Musashi homolog 1.
Human Swiss Prot No.:	O43347
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Musashi-1 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year