

Antibodies Product Datasheet

Product Name:	Cytochrome P450 17A1 mouse mAb
Catalog Number:	N0904
Reactivity:	Human
Applications:	WB
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human Cytochrome P450 17A1 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	cyp17a1
Protein Name:	
Other Name:	20 lyase;CP17A_HUMAN;CPT7;CYP17;CYP17A1;CYPXVII;Cytochrome P450 17A1;Cytochrome P450 family 17;Cytochrome P450 family 17 subfamily A polypeptide 1;Cytochrome p450 subfamily XVII (steroid 17 alpha hydroxylase) adrenal hyperplasia;Cytochrome p450 XVIIA1; Cy
Human Swiss Prot No.:	P05093
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Cytochrome P450 17A1 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year