

Antibodies Product Datasheet

Product Name:	CrkII mouse mAb
Catalog Number:	N0906
Reactivity:	Human
Applications:	WB
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human CrkII protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	crk
Protein Name:	
Other Name:	Adapter molecule crk;Avian sarcoma virus CT10 (v crk) oncogene homolog;CRK;CRK isoform 2;CRK isoform II;CRK_HUMAN;CRKII;FLJ38130;OTTHUMP00000115366;OTTHUMP0000198330;p38; Proto oncogene C crk;Proto-oncogene C-cr;v crk avian sarcoma virus CT10 oncogene h
Human Swiss Prot No.:	P46108
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CrkII and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year