

Antibodies Product Datasheet

Product Name:	CBR1 mouse mAb
Catalog Number:	N0907
Reactivity:	Human
Applications:	WB, ICC
Dilutions:	wb 1:1000 icc 1:100
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human CBR1 protein fragments expressed in E.coli
Concentration:	1 mg/ml
Gene Name:	cbr1
Protein Name:	
Other Name:	15 hydroxyprostaglandin dehydrogenase [NADP]; 15-hydroxyprostaglandin dehydrogenase [NADP]; Carbonyl reductase [NADPH] 1; CBR 1; CBR1; CBR1_HUMAN; CRN; NADPH dependent carbonyl reductase 1; NADPH-dependent carbonyl reductase 1; Prostaglandin 9 ketoreductase; Pro
Human Swiss Prot No.:	P16152
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CBR1 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year