

Antibodies Product Datasheet

Product Name:	clAP2 mouse mAb
Catalog Number:	N0908
Reactivity:	Human, Monkey
Applications:	WB
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human c-IAP2 protein fragments expressed in E.coli
Concentration:	1 mg/ml
Gene Name:	birc3
Protein Name:	
Other Name:	AIP 1;AIP1;API 2;API2;API2;Apoptosis inhibitor 2;Baculoviral IAP repeat containing 3;Baculoviral IAP repeat containing protein 3;Baculoviral IAP repeat-containing protein 3;BIRC 3;BIRC3;BIRC3;BIRC3_HUMAN;C IAP2;C-IAP2;CIAP 2;CIAP 2;CIAP2;HAIP 1;HAIP1;HAIP
Human Swiss Prot No.:	Q13489
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of c-IAP2 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year