

Antibodies Product Datasheet

Product Name:	Placental alkaline phosphatase (PLAP) mouse mAb
Catalog Number:	N0914
Reactivity:	Human
Applications:	WB
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant Placental alkaline phosphatase (PLAP) protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	alpp
Protein Name:	
Other Name:	Alkaline phosphatase; Alkaline phosphatase placental; Alkaline phosphatase placental type; Alkaline phosphatase Regan isozyme; ALP; Alp1; ALPP; FLJ61142; Germ-cell alkaline phosphatase; nagao isozyme; OTTHUMP00000164354; PALP; Placental alkaline phosphata
Human Swiss Prot No.:	P05187
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Placental alkaline phosphatase (PLAP) and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year