

Antibodies Product Datasheet

Product Name:	CHD4 mouse mAb
Catalog Number:	N0918
Reactivity:	Human
Applications:	WB,IP
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human CHD4 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	chd4
Protein Name:	
Other Name:	ATP dependent helicase CHD4;ATP-dependent helicase CHD4;CHD 4;CHD-4;CHD4;CHD4_HUMAN; Chromodomain helicase DNA binding protein 4;Chromodomain-helicase-DNA-binding protein 4;DKFZp686E06161;Mi 2 autoantigen 218 kDa protein;Mi 2b;Mi-2 autoantigen 218 kDa pro
Human Swiss Prot No.:	Q14839
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CHD4 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year