

Antibodies Product Datasheet

Product Name:	LIN28A mouse mAb
Catalog Number:	N0920
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human LIN28A protein fragments expressed in E.coli
Concentration:	1 mg/ml
Gene Name:	lin28a
Protein Name:	
Other Name:	CSDD1; FLJ12457; LIN 28; Lin-28A; LIN28; LIN28A; LN28A_HUMAN; Protein lin-28 homolog A; ZCCHC1; Zinc finger CCHC domain-containing protein 1.
Human Swiss Prot No.:	Q9H9Z2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of LIN28A and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year