

Antibodies Product Datasheet

Product Name:	Annexin VI mouse mAb
Catalog Number:	N0921
Reactivity:	Human
Applications:	WB
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human Annexin VI protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	anxa6
Protein Name:	
Other Name:	67 kDa calelectrin;Annexin A6;Annexin VI;Annexin VI p68;Annexin-6;AnnexinA6;AnnexinVI; ANX 6;ANX A6;ANX6;ANXA6;ANXA6_HUMAN;Calcium binding protein p68;Calelectrin;Calphobindin II;Calphobindin-II;CalphobindinII;CBP 68;Chromobindin 20;CPBII;Lipocortin VI; L
Human Swiss Prot No.:	P08133
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Annexin VI and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year