

Antibodies Product Datasheet

Product Name:	DR5 mouse mAb
Catalog Number:	N0923
Reactivity:	Human, Mouse
Applications:	WB, ICC
Dilutions:	wb 1:500-1:2000 icc 1:100
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human DR5 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	tnfrsf10b
Protein Name:	
Other Name:	Fas like protein; Apoptosis inducing protein TRICK2A/2B; Apoptosis inducing receptor TRAIL R2; CD 262; CD262; CD262 antigen; Cytotoxic TRAIL receptor 2; Death domain containing receptor for TRAIL/Apo 2L; Death domain containing receptor for TRAIL/Apo2L; Death rece
Human Swiss Prot No.:	O14763
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DR5 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year