

Antibodies Product Datasheet

Product Name:	IL18 binding protein mouse mAb
Catalog Number:	N0929
Reactivity:	Human
Applications:	WB
Dilutions:	wb 1:400
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human IL18 binding protein protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	il18bp
Protein Name:	
Other Name:	I18BP_HUMAN;IL-18BP;IL18 BP;IL18 BPα;IL18BP;IL18BPα;Interleukin 18 binding protein; Interleukin-18-binding protein;MC51L 53L 54L homolog gene product;Tadekinig alfa; Tadekinig-alfa.
Human Swiss Prot No.:	O95998
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of IL18 binding protein and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year