

Antibodies Product Datasheet

Product Name:	ALDH2 mouse mAb
Catalog Number:	N0931
Reactivity:	Human, Mouse, Rat, Monkey, Hamster
Applications:	WB
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human ALDH2 protein fragments expressed in E.coli
Concentration:	1 mg/ml
Gene Name:	aldh2
Protein Name:	
Other Name:	Acetaldehyde dehydrogenase 2; Aldehyde dehydrogenase 2 family (mitochondrial); Aldehyde dehydrogenase 2 family; Aldehyde dehydrogenase mitochondrial; Aldehyde dehydrogenase, mitochondrial; ALDH 2; ALDH class 2; ALDH E2; ALDH-E2; Aldh2; ALDH2_HUMAN; ALDHI; ALDM; Liver
Human Swiss Prot No.:	P05091
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ALDH2 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year