

Antibodies Product Datasheet

Product Name:	DDB1 mouse mAb
Catalog Number:	N0932
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human DDB1 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	ddb1
Protein Name:	
Other Name:	Damage specific DNA binding protein 1; Damage-specific DNA-binding protein 1; DDB 1; DDB p127 subunit; Ddb1; DDB1_HUMAN; DDBa ; DNA damage binding protein 1; DNA damage-binding protein 1; DNA damage-binding protein a; HBV X-associated protein 1; UV damaged DNA bindi
Human Swiss Prot No.:	Q16531
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DDB1 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year