

Antibodies Product Datasheet

Product Name:	IMP3 mouse mAb
Catalog Number:	N0933
Reactivity:	Human
Applications:	WB
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human IMP3 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	igf2bp3
Protein Name:	
Other Name:	Cancer/testis antigen 98;CT98;DKFZp686F1078;hKOC;IF2B3_HUMAN;IGF II mRNA binding protein 3;IGF-II mRNA-binding protein 3;IGF2 mRNA binding protein 3;IGF2 mRNA-binding protein 3;IGF2BP3;IMP 3;IMP-3;Insulin like growth factor 2 mRNA binding protein 3;Insuli
Human Swiss Prot No.:	O00425
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of IMP3 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year