

Antibodies Product Datasheet

Product Name:	FAK mouse mAb
Catalog Number:	N0934
Reactivity:	Human
Applications:	WB,IF/ICC
Dilutions:	wb 1:1000 icc 1:200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human FAK protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	ptk2
Protein Name:	
Other Name:	FADK 1;FADK;FAK 1;FAK 1;FAK related non kinase polypeptide;FAK1;FAK1_HUMAN;Focal adhesion kinase 1;Focal adhesion kinase 1;Focal adhesion Kinase;Focal adhesion kinase isoform FAK Del33;Focal adhesion kinase related nonkinase;FRNK;p125FAK;pp125FAK;pp125FAK
Human Swiss Prot No.:	Q05397
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FAK and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year