

Antibodies Product Datasheet

Product Name:	Cystatin C mouse mAb
Catalog Number:	N0936
Reactivity:	Human
Applications:	Elisa
Dilutions:	ELISA 1:10000-20000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant human cystatin c protein.
Concentration:	1 mg/ml
Gene Name:	cystatin c
Protein Name:	
Other Name:	AD 8;AD8;Amyloid angiopathy and cerebral hemorrhage;ARMD11;bA218C14.4 (cystatin C);bA218C14.4;Cst 3;Cst3;CST3 protein;Cystatin 3;Cystatin-3;Cystatin-C;Cystatin3;CystatinC;CYTC_HUMAN;Gamma trace;Gamma-trace;HCCAA;MGC117328;Neuroendocrine basic polypeptide;
Human Swiss Prot No.:	P01034
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	
Storage and Stability:	-20°C/1 year