

Antibodies Product Datasheet

Product Name:	Cystatin C mouse mAb
Catalog Number:	N0939
Reactivity:	Human
Applications:	sELISA, Detector
Dilutions:	ELISA 1:10000-20000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant human cystatin c protein.
Concentration:	1 mg/ml
Gene Name:	cyc
Protein Name:	
Other Name:	CST3;Cystatin 3.
Human Swiss Prot No.:	P01034
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	
Storage and Stability:	-20°C/1 year