

Antibodies Product Datasheet

Product Name:	CD31 mouse mAb
Catalog Number:	N0940
Reactivity:	Human
Applications:	WB,IP
Dilutions:	wb dilution 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human CD31 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	pecam1
Protein Name:	
Other Name:	Adhesion molecule;CD31;CD31 antigen;CD31 EndoCAM;Endocam;FLJ34100;FLJ58394;GPIIA; GPIIA';PECA1;PECA1_HUMAN;Pecam 1;PECAM 1 CD31 EndoCAM;PECAM;PECAM-1;Pecam1;Platelet endothelial cell adhesion molecule;Platelet/endothelial cell adhesion molecule 1; Platele
Human Swiss Prot No.:	P16284
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CD31 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year