

Antibodies Product Datasheet

Product Name:	WDR77 mouse mAb
Catalog Number:	N0949
Reactivity:	Human, Mouse, Rat
Applications:	WB, ICC
Dilutions:	wb dilution 1:1000 icc dilution 1:100
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human WDR77 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	mep50/wdr77
Protein Name:	
Other Name:	2610312E17Rik; Androgen receptor cofactor p44; C79984; HKMT1069; MEP 50; MEP-50; MEP50; MEP50_HUMAN; Methylosome protein 50; MGC2722; Nbla10071; p44; p44/Mep50; RGD1310479; RP11 552M11.3; WD repeat containing protein 77; WD repeat domain 77; WD repeat-containing protein
Human Swiss Prot No.:	Q9BQA1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of WDR77 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year