

Antibodies Product Datasheet

Product Name:	Bax mouse mAb
Catalog Number:	N0952
Reactivity:	Human, Monkey
Applications:	WB
Dilutions:	wb dilution 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human Bax protein fragments expressed in E.coli
Concentration:	1 mg/ml
Gene Name:	bax
Protein Name:	
Other Name:	Apoptosis regulator BAX; Apoptosis regulator BAX cytoplasmic isoform beta; Apoptosis regulator BAX membrane isoform alpha; BAX; Bax isoform psi; BAX protein cytoplasmic isoform delta; Bax protein cytoplasmic isoform delta.; Bax protein cytoplasmic isoform
Human Swiss Prot No.:	Q07812
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Bax and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year