

Antibodies Product Datasheet

Product Name:	RPA32/RPA2 mouse mAb
Catalog Number:	N0953
Reactivity:	Human
Applications:	WB,IP,ICC
Dilutions:	wb dilution 1:2000 icc dilution 1:200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human RPA32/RPA2 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	rpa2
Protein Name:	
Other Name:	60S acidic ribosomal protein P1;AA409079;AI325195;AU020965;HSSB;ik:tdsubc_2g1;M(2)21C; MGC137236;OTTHUMP00000004008;p32;p34;RCJMB04_6d17 replication protein A2, 32kDa;REPA 2; REPA1;REPA2;Replication factor A protein 2;Replication protein A 32 kDa subunit;
Human Swiss Prot No.:	P15927
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of RPA32/RPA2 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year