

Antibodies Product Datasheet

Product Name:	BRAP mouse mAb
Catalog Number:	N0958
Reactivity:	Human
Applications:	WB
Dilutions:	wb dilution 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human BRAP protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	brap
Protein Name:	
Other Name:	3010002G07Rik;BRAP2;BRCA1 associated protein;EC 6.3.2.;Galectin 2 binding protein;IMP; Impedes mitogenic signal propagation;Renal carcinoma antigen NY REN 63; RING finger protein 52;RNF52;zgc:92894.
Human Swiss Prot No.:	Q7Z569
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of BRAP and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year