

Antibodies Product Datasheet

Product Name:	hHR23b mouse mAb
Catalog Number:	N0959
Reactivity:	Human, Mouse, Rat, Monkey, Hamster
Applications:	WB, ICC, IHC
Dilutions:	wb dilution 1:1000 icc dilution 1:100
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human hHR23b protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	hr23b
Protein Name:	
Other Name:	hHR 23b; hHR23B; HR 23B; HR23B; mHR 23B; mHR23B; P58; RAD 23B; RAD23 (S. cerevisiae) homolog B; RAD23 homolog B (S. cerevisiae); RAD23 homolog B; RAD23 yeast homolog of B; RAD23B; RD23B_HUMAN; UV excision repair protein RAD23 homolog B; XP C repair complementing compl
Human Swiss Prot No.:	P54727
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of hHR23b and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year