

## Antibodies Product Datasheet

|                               |                                                                                                                                                                                                                                                                 |
|-------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | HTRA2 mouse mAb                                                                                                                                                                                                                                                 |
| <b>Catalog Number:</b>        | N0960                                                                                                                                                                                                                                                           |
| <b>Reactivity:</b>            | Human                                                                                                                                                                                                                                                           |
| <b>Applications:</b>          | WB                                                                                                                                                                                                                                                              |
| <b>Dilutions:</b>             | wb dilution 1:500                                                                                                                                                                                                                                               |
| <b>Modification:</b>          | Unmodified/ Total                                                                                                                                                                                                                                               |
| <b>Source:</b>                | Monoclonal Mouse                                                                                                                                                                                                                                                |
| <b>Purification:</b>          | The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.                                                                                                                                              |
| <b>Immunogen:</b>             | Recombinant human HtrA2/Omi protein.                                                                                                                                                                                                                            |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                                                                                                                                                                         |
| <b>Gene Name:</b>             | HtrA2/Omi                                                                                                                                                                                                                                                       |
| <b>Protein Name:</b>          |                                                                                                                                                                                                                                                                 |
| <b>Other Name:</b>            | High temperature requirement protein A2;HTRA 2;HtrA like serine protease;HtrA serine peptidase 2;HtrA, E. coli, homolog of, 2;HtrA2;HTRA2_HUMAN;mitochondrial;Omi stress regulated endoprotease;Omi stress-regulated endoprotease;PARK 13;PARK13;Protease serin |
| <b>Human Swiss Prot No.:</b>  | O43464                                                                                                                                                                                                                                                          |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                                                                                                         |
| <b>Specificity:</b>           | This antibody detects endogenous levels of HtrA2/Omi and does not cross-react with related proteins.                                                                                                                                                            |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                                                                                                                                                                    |