

Antibodies Product Datasheet

Product Name:	Hsp90 mouse mAb
Catalog Number:	N0963
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	wb dilution 1:2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human HSP90AA1 protein fragments expressed in E.coli
Concentration:	1 mg/ml
Gene Name:	hsp90
Protein Name:	
Other Name:	D6S182; FLJ26984; FLJ31884; Heat shock 86 kDa; heat shock 90kDa protein 1 alpha; Heat shock protein 90kDa alpha cytosolic class A member 1; Heat shock protein 90kDa alpha cytosolic class B member 1; Heat shock protein HSP 90 alpha; Heat shock protein HSP 90 beta
Human Swiss Prot No.:	P07900
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Hsp90 alpha and Hsp90 beta.
Storage and Stability:	-20°C/1 year