

Antibodies Product Datasheet

Product Name:	Hsp90 β mouse mAb
Catalog Number:	N0964
Reactivity:	Human, Mouse, Rat, Monkey, Hamster
Applications:	WB, ICC, IP
Dilutions:	wb dilution 1:2000 icc dilution 1:200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human Hsp90 beta protein fragments expressed in E.coli
Concentration:	1 mg/ml
Gene Name:	hsp90ab1
Protein Name:	
Other Name:	HSP90AB1; 90 kda heat shock protein beta HSP90 beta; D6S182; FLJ26984; Heat shock 84 kDa; Heat shock 90kD protein 1, beta; Heat shock 90kDa protein 1 beta; Heat shock protein 90kDa alpha (cytosolic) class B member 1; Heat shock protein beta; Heat shock protein HSP
Human Swiss Prot No.:	P08238
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Hsp90 beta and does not cross-react with Hsp90 alpha and other proteins.
Storage and Stability:	-20°C/1 year