

Antibodies Product Datasheet

Product Name:	Ku70 mouse mAb
Catalog Number:	N0965
Reactivity:	Human, Monkey
Applications:	WB, ICC, IP
Dilutions:	wb dilution 1:1000 icc dilution 1:200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human Ku70 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	xrcc6
Protein Name:	
Other Name:	5"-deoxyribose-5-phosphate lyase Ku70; 5"-dRP lyase Ku70; 70 kDa subunit of Ku antigen; ATP dependent DNA helicase 2 subunit 1; ATP dependent DNA helicase II 70 kDa subunit; ATP-dependent DNA helicase 2 subunit 1; ATP-dependent DNA helicase II 70 kDa subunit
Human Swiss Prot No.:	P12956
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Ku70 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year