

Antibodies Product Datasheet

Product Name:	JNK1 mouse mAb
Catalog Number:	N0968
Reactivity:	Human, Mouse, Rat
Applications:	WB, IF/ICC
Dilutions:	wb dilution 1:1000 icc dilution 1:100
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human JNK1 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	mapk8
Protein Name:	
Other Name:	AI849689 ; c Jun N terminal kinase 1 ; C-JUN kinase 1 ; c-Jun N-terminal kinase 1 ; EC 2.7.11.24 ; JAK 1A ; JAK1A ; JNK 1 ; JNK 46 ; JNK ; JNK-46 ; JNK1A2 ; JNK21B1/2 ; MAP kinase 8 ; MAPK 8 ; MAPK8 ; Mitogen activated protein kinase 8 ; Mitogen-activated protein kinase 8
Human Swiss Prot No.:	P45983
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of JNK1 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year