

Antibodies Product Datasheet

Product Name:	Ku80 mouse mAb
Catalog Number:	N0969
Reactivity:	Human, Monkey
Applications:	WB, IF/ICC, IP
Dilutions:	wb dilution 1:1000 icc dilution 1:400 ip dilution 1:100
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human Ku80 protein fragments expressed in E.coli
Concentration:	1 mg/ml
Gene Name:	xrcc5
Protein Name:	
Other Name:	ATP dependent DNA helicase II 80 kDa subunit; ATP dependent DNA helicase II 86 Kd subunit; ATP dependent DNA helicase II; ATP-dependent DNA helicase 2 subunit 2; ATP-dependent DNA helicase II 80 kDa subunit; CTC box binding factor 85 kDa; CTC box-binding
Human Swiss Prot No.:	P13010
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Ku80 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year