

Antibodies Product Datasheet

Product Name:	EGF Receptor mouse mAb
Catalog Number:	N0972
Reactivity:	Human, Monkey
Applications:	WB, ICC, IP
Dilutions:	wb dilution 1:1000 icc dilution 1:200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human EGFR protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	egfr
Protein Name:	
Other Name:	Avian erythroblastic leukemia viral (v erb b) oncogene homolog; Cell growth inhibiting protein 40; Cell proliferation inducing protein 61; EGF R; EGFR; EGFR_HUMAN; Epidermal growth factor receptor (avian erythroblastic leukemia viral (v erb b) oncogene homolog)
Human Swiss Prot No.:	P00533
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous level of total EGFR and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year